Montage- und Gebrauchsanweisung	04	DEU
Installation and operating instructions	15	ENG
Mode d'emploi	26	FRA
Istruzioni per il montaggio e l'uso	37	ITA
Instrucciones de montaje y funcionamiento	48	ESP
Montage- en gebruiksaanwijzing	59	NLD
Bruksanvisning	70	SWE
取り付けおよび使用の説明書	81	JAP

- D^{med}® **OPTICLUX** 10-1P TX
- D^{med}® **OPTICLUX** 10-1 C T1
- D^{med}® **OPTICLUX** 10-2 P TX
- D^{med}® **OPTICLUX** 10-2 C T1

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Examination light

Lampe d'examens

Apparecchio per visita

Luz de examen

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Undersökningsarmatur

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SYMBOLE; SYMBOLS; SYMBOLES; SIMBOLI; SÍMBOLOS; SYMBOLEN; SYMBOLER; 記号の意味

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IMPORTANT!
THIS OPERATING MANUAL MUST BE THOROUGHLY READ BEFORE
USING THE PRODUCT!


→ **STORE FOR LATER USE!**

CONTENT

1.	VARIANTS AND SCOPE OF DELIVERY	16
1.1	OPTICLUX 10-1 C T1 / OPTICLUX 10-2 C T1	16
1.2	OPTICLUX 10-1 P TX / OPTICLUX 10-2 P TX.....	16
2.	SAFETY INSTRUCTIONS	17
2.1	Intended use	17
2.2	User profiles	17
2.3	Safety instructions	17
2.4	Warning levels	17
3.	ASSEMBLY of OPTICLUX 10-1 C T1 OPTICLUX 10-2 C T1	18
3.1	Load data.....	18
3.2	Shortening the ceiling tube	18
3.3.	Installing the ceiling holder	18
3.4.	Installing the ceiling tube	19
3.5	Assemble the light on ceiling tube	19
4.	ASSEMBLY OPTICLUX 10-1 P TX OPTICLUX 10-2 P TX	20
4.1	Load data.....	20
4.2	Assembly	20
5.	OPERATION.....	20
5.1	OPTICLUX 10-1 C T1 / OPTICLUX 10-1 P TX.....	20
5.2	OPTICLUX 10-2 C T1 / OPTICLUX 10-2 P TX.....	20
5.3	Transport position	20
6.	CLEANING AND DISINFECTION.....	21
7.	SAFETY INSPECTIONS	21
8.	DISMANTLING	21
8.1	Disposal.....	21
9.	ACCESSORIES.....	22
10.	ADDITIONAL NOTES.....	22
11	TROUBLESHOOTING.....	23
12	TECHNICAL DATA.....	23
13.	ELECTROMAGNETIC COMPATIBILITY (EMC)	24

1. VARIANTS AND SCOPE OF DELIVERY

1.1 OPTICLUX 10-1 C T1 / OPTICLUX 10-2 C T1

A: Ceiling plate and canopy


B: Retaining ring

C: Ceiling tube

E: Magnifier cover

D: Arm system with lamp head

1.2 OPTICLUX 10-1 P TX / OPTICLUX 10-2 P TX

A: Arm system with lamp head

B: Magnifier cover

C: Power cable

2. SAFETY INSTRUCTIONS

2.1 Intended use

The light Dmed® OPTICLUX is an examination light. It is intended to locally illuminate a patient's body to support dermatological diagnosis and treatment, at close range of up to 15 cm from the surface to be treated. An interruption of the diagnosis or treatment caused by a light failure is always possible without posing any risk for the patient. The light cannot be used in operating rooms.

Moreover, Dmed® OPTICLUX can also be used in the medical environment of laboratory or in quality control.

2.2 User profiles

Health Professional

All individuals who have completed medical training and work in the professional field they trained for.

Cleaning specialist

Trained in national and workplace hygiene regulations.

Qualified electrician

Trained in electronics and electrical technology and knows the relevant standards and regulations.

Qualified specialist

Qualified due to his technical training, knowledge and experience and knowledge of the rules, to perform assembly/dismantling.

2.3 Safety instructions

- ▶ Operation by health professionals
- ▶ This manual is part of the product and must be stored and made available to all future users.
- ▶ All work on the luminaire (including repairs) may only be performed by a qualified electrician. Assembly may only be done by a qualified specialist.
- ▶ The luminaire may not be altered or tampered with. Only approved original parts may be used. Any use other than that intended, using original parts, can change technical parameters and pose a death hazard.
- ▶ Do not exceed the maximum weight and do not hang, lean or climb on the arm, because this may cause the device to tip over, which could result in serious injury.
- ▶ Operation in potentially explosive areas is prohibited. The luminaire's current supply is a potential ignition source.
- ▶ The luminaire must be used only in dry, dust-free rooms.
- ▶ The luminaire should not be left on without supervision.
- ▶ Only connect the luminaire to the supply network with a protective ground conductor (PE), to prevent electrical shock.
- ▶ In the case of luminaires of protection class I, the protective ground conductor (PE) must be connected to luminaires housing.
- ▶ Do not use a damaged luminaire. Defective cables and a defective hand grip also pose a potential hazard. Do not place a cord near heat sources or sharp edges.

- ▶ Do not place extra loads on the luminaire head and arm system.
- ▶ The luminaire must not be covered with a cloth or similar item during operation.
- ▶ The ventilation openings (if present) must always be kept clear during operation!
- ▶ The luminaire must not be operated near external heat sources that exceed the luminaire's maximum ambient temperature.
- ▶ The luminaire must not be used outside the specified ambient conditions.
- ▶ Do not use with medical devices that may react sensitively to a light spectrum within the visible range (such as pulsating light and/or light with high illumination intensity).
- ▶ The luminaire may only be used for the intended use described here.
- ▶ The manufacturer cannot be held responsible for any damages resulting from use deviating from the intended use, or failure to observe the safety instructions and warnings.
- ▶ When using more than one luminaire at the same time, the total illumination intensity must not exceed $E_e 1000W/m^2$.
- ▶ Before connecting to the power grid, it is essential to check that the grid data match the device data.
- ▶ The luminaire must be secured during transport within the clinic.

2.4 Warning levels

DANGER

Warning of hazards that can result in **death or serious injury** if instructions are not followed.

WARNING

Warning of hazards that can result in **injury** if instructions are not followed.

CAUTION


Warning of hazards that can cause **material damage** if instructions are not followed.

3. ASSEMBLY of OPTICLUX 10-1 C T1 OPTICLUX 10-2 C T1


3.1 Load data

Bending moment M_B	25Nm
Vertical weight F_G	90N

3.2 Shortening the ceiling tube

- ▶ Before shortening the ceiling tube, remove the cable from the ceiling tube.
- ▶ Use a hack saw to shorten the upper end of the ceiling tube to the desired length, and then deburr it.

- ▶ Remove the fastening screw "C".
- ▶ Insert the ceiling tube "A" into the ceiling bracket "B" and drill the ceiling holder's existing hole to 9 mm in diameter. Drill the opposite hole separately.
- ▶ Note: Pull the cable through from the lower side to the upper pipe side (the 3-pin plug first) after sawing and drilling.

3.3. Installing the ceiling holder

- ▶ **Fastening hardware** does **not** come with delivery.

⚠ DANGER

Assembly by qualified personnel

- ▶ Assembly must be done by qualified personnel only. Lack of appropriate knowledge may pose a death hazard.
- ▶ Two people are needed for assembly.

⚠ DANGER


Death hazard from falling luminaire.

- ▶ The ceiling mount must only be mounted on ceilings with a concrete stability class of B25 (C20/25) or higher.
- ▶ There should be no contact with reinforcement parts of the solid ceiling during mounting. When in doubt, a licensed technician must approve installation on the specific installation surface. The load bearing capacity of the ceiling structure must be planned, checked and confirmed by a structural engineer beforehand.
- ▶ The holes must be drilled by a professional and in compliance with the drilling tolerances permitted by the manufacturer of the fastening anchor. If a drilling error occurs, such as drilling into a reinforcement bar, a structural engineer must be consulted.
- ▶ Install the luminaire such that the vertical stops are not continuously stressed during operation.
- ▶ With plaster or cladding in front of concrete, the fastening anchor must be driven completely into the concrete.
- ▶ The screws must be tightened carefully using a torque wrench according to the specifications of the fastening anchor manufacturer.

⚠ DANGER

Death hazard from electric shock.


- ▶ All poles of the power cable must be disconnected from the power grid from an external lockable switch.


- ▶ Draw 6 drill marks.
- ▶ Observe the position of opening "A" for connecting the current.

WARNING


Wear safety equipment according to the tool manufacturer's instructions.

- ▶ Drill holes and blow out with a bellows.

- ▶ Check the hole distances.


- ▶ Hold the ceiling bracket to the ceiling and hammer in the fastening anchor.
- ▶ Tighten the fastener according to manufacturer's instructions.

DANGER


Death hazard from electric shock.

- ▶ Do not turn on the power supply before the luminaire head has been installed.
- ▶ This device must only be connected to a power grid with a protective ground conductor (PE) to prevent risk of electric shock


- ▶ Connect the power.

3.4. Installing the ceiling tube


- ▶ Pull the ceiling tube cable through the opening "B" of the ceiling holder.
- ▶ Insert the ceiling tube "A" into the ceiling holder.
- ▶ Secure the safety screw "C" and the M8 nut.

- ▶ Tighten the M8 safety screw "C" and the nut (20 Nm).
- ▶ Tighten all 4 grub screws "D" (5 Nm).

- ▶ Connect the secondary cable to the connection point
- ▶ Pay attention to the polarity:
+ = black
- = white


- ▶ Slide the ceiling cover and ring over the ceiling bracket and screw them tight.

3.5 Assemble the light on ceiling tube

WARNING

The light will fall if fastening is done incorrectly

- ▶ The securing element must be screwed in fully and correctly.

- ▶ Insert the light and the adapter 1 into the ceiling tube 2
- ▶ Insert securing element with threaded pin 3 and screw in to the stop
- ▶ Loosen by one quarter-turn to ensure the operation of the turnover axis

4. ASSEMBLY OPTICLUX 10-1 P TX OPTICLUX 10-2 P TX

4.1 Load data

Bending moment M_B	25Nm
Vertical weight F_G	20N

4.2 Assembly

- ▶ The lights are equipped with an adapter pin. The light must be positioned in one of the accessories described in chapter 9.

5. OPERATION

CAUTION


- ▶ This product emits potentially hazardous optical radiation. Do not stare directly into the light cone. Eye irritation may occur.
- ▶ The radiation emitted by this product complies with the exposure limits for reducing the risk of photobiological hazards based on IEC 62471: RG 1 (low risk).

DANGER

Death hazard from electric shock.


- ▶ Do not plug in any damaged power cables.
- ▶ If the power cable shows damage, immediately replace it with a new one.
- ▶ The connection voltage and frequency must match the data on the type plate.
- ▶ Only connect to power grid with protective ground conductor (PE)

OPTICLUX 10-1 C T1 / OPTICLUX 10-2 C T1


- ▶ Insert power cable
- ▶ Connect the cable to the mains
- ▶ Perform a function test before every use:
- ▶ All LED's in the light cone must illuminate.

5.1 OPTICLUX 10-1 C T1 / OPTICLUX 10-1 P TX

- ▶ Switch the LED module on/off
- ▶ Press and hold the button for dimming

5.2 OPTICLUX 10-2 C T1 / OPTICLUX 10-2 P TX

- ▶ A: Switch the LED module on/off
- ▶ A: Press and hold the button for dimming
- ▶ B: Switch to Woodlight

5.3 Transport position

- ▶ Only transport the light in this position.

6. CLEANING AND DISINFECTION

DANGER

Death hazard from electric shock

- ▶ Before disinfection cleaning, switch off the power connection and secure against unintentional activation.

CAUTION

Material damage due to incorrect cleaning

- ▶ For cleaning, use only products that do not impair operation of the luminaire.
- ▶ For cleaning, do not use solvent-based, chlorine-based or abrasive detergents, because they can crack plastic parts and cause other damage.
- ▶ The cleaning agents must be approved for use on plastics such as PC, PMMA, PA and ABS.
- ▶ Concentrated disinfectant can damage the luminaire.
- ▶ For concentration and application times, check the information provided with the product used.
- ▶ The wrong cloth may cause scratches.

RECOMMENDED DISINFECTANTS

- | | |
|-----------------------------|--------------------|
| ▶ Bacillol 30 Foam | ▶ Terralin PAA |
| ▶ Dismozon Plus | ▶ Terralin Protect |
| ▶ Kohrsolin Extra | ▶ Virex Tb |
| ▶ Lysoformin | ▶ CaviCide 1 |
| ▶ Mediclean | ▶ Gemicidal Bleach |
| ▶ Meliseptol Foam Pure | ▶ Hexaquart XL |
| ▶ Microbac Tissues | |
| ▶ Mikrozyd Sensitive Liquid | |

CAUTION

Dirt reduces the light strength

- ▶ Keep cover clear through regular cleaning.
- ▶ Only wipe cleaning allowed.

- ▶ Clean the PMMA screen with a non-abrasive cleaning cloth and a suitable cleaning agent.

CAUTION

To minimize the risk of disease transmission, applicable health and safety regulations and the requirements of the national hygiene and disinfection authorities must be observed in addition to these instructions.

7. SAFETY INSPECTIONS

DANGER

Death hazard from electric shock.

- ▶ Unplug from the grid.
- ▶ The power supply cable must be checked at least once a year for damage.

CAUTION

- ▶ Maintenance and repairs can be performed only by qualified electricians.
- ▶ The corresponding user profile is in Chapter 2 Safety instructions.

EVERY YEAR:

- ▶ Check the power supply cable for damage and replace if necessary.
- ▶ Check for deformations and cracks in the plastic parts.
- ▶ Check the load bearing system for distortion or damage.
- ▶ Check for loose parts.

8. DISMANTLING

DANGER

Death hazard from electric shock.

Before dismantling the luminaire, disconnect all pins from the power supply.

8.1 Disposal

Do not discard the luminaire with household waste. Follow local regulations and bring the luminaire to a disposal site or give it to a dealer with appropriate service.
Cut the cable directly at the casing.

The products listed above are over 95% recyclable. For a high percentage of the used materials to be either physically reused or used for energy after the end of this product's life cycle, the luminaires have been designed with recycling in mind. They do not contain materials that are hazardous or require monitoring.

9. ACCESSORIES

Table clamp
(Art. No.: D14.228.000 - black)

Universal mount
(Art. No.: D13.430.000 - pure white)

Rail clamp
(Art. No.: D13.269.000)

Wall bracket
(Art. No.: D13.231.000 - pure white)

Roller stand
(Art. No.: D15.595.000)

10. ADDITIONAL NOTES

The luminaire itself is maintenance free.

Additional documents may be requested from the manufacturer for this product.

Using this luminaire does not pose a risk to other equipment.

To save energy, the luminaire should be switched on only when it is actually needed.

Any serious incident that has occurred with the product **must be reported** to the manufacturer or their representative and the responsible authorities of the country in which the user is located.

11 TROUBLESHOOTING

Fault	Possible cause	Troubleshooting	User profiles
The light does not light	Contact problem	Switch on again	All
The light does not light	No mains voltage	Check mains voltage, /all connections	Electrician
The light does not light	LED module defective	Contact manufacturer service	By manufacturer service only

12 TECHNICAL DATA

Electrical data:		
Rated input voltage	100 – 240 V AC	
Frequency range	50-60 Hz	
Power consumption		
OPTICLUX 10-1 C T1 / OPTICLUX 10-1 P TX	5-10 VA	
OPTICLUX 10-2 C T1 / OPTICLUX 10-2 P TX	11-17 VA	
Input current	0.03-0.11 A	
Power factor	0.52-0.97	
Power supply secondary side	24VDC	
Photometric values*:		
Central illuminance E_v at 15 cm (0.5 feet) distance (380 – 780nm, 6500K)	6'000 lx	
Central irradiance E_e at 15 cm (0.5 feet) distance (315 – 400nm, Wood light)	7,6 W/m ²	
Light field diameter d10 at 15 cm (0.5 feet) distance	Ø = 56 cm (1.9")	
Colour temperature	6500 K	
Colour rendering index Ra	>93	
Colour rendering index R9	>90	
Total irradiance E_e at maximum intensity	<25 W/m ²	
Ratio irradiance E_e to illuminance E_c	4 mW/m ² /lx	* -10% / +20% tolerance
Ambient conditions for transport, storage and operation:		
Ambient temperature (storage and transport)	-20°C to +70°C (-4°F to +158°F)	
Ambient temperature (operation)	+10°C to +35°C (+50°F to +95°F)	
Relative humidity (non-condensing) (storage and transport)	max. 90%	
Relative humidity (non-condensing) (operation)	max. 75%	
Weight:		
OPTICLUX 10-1 C T1 / OPTICLUX 10-2 C T1	9kg (19.8 lb)	
OPTICLUX 10-1 P TX / OPTICLUX 10-2 P TX	2kg (4.4 lb)	
Mode of operation:		
Mode of operation	Continuous operation	

Classification:	
OPTICLUX 10-1 C T1 / OPTICLUX 10-2 C T1 OPTICLUX 10-1 P TX / OPTICLUX 10-2 P TX	Protection class I (C-version) Protection class II (P-version)
Degree of protection according to IEC 60529	IP 20
Classification according to EU REGULATION 2017/745 (MDR), article 51 U.S. FDA Device Class	Class I Class I
Electrical safety testing and EMC according to:	AAMI ES60601-1 : 2005/A2 :2010/(R)2012 CAN/CSA-C22.2 No 60601-1:14 IEC 60601-1:2005 + A1 :2012 IEC 60601-1-2:2014 IEC 60601-2-41:200- + A1 :2013
Blue light hazard according to EN/IEC 62471	RG 1 (low risk)
Service life of the light source:	
Service life LED	50'000h L80/B10
Service life UV-LED	30'000h L70/B50

13. ELECTROMAGNETIC COMPATIBILITY (EMC)

Electrical medical devices are subject to special precautionary measures regarding electromagnetic compatibility. This device can be affected by other electrical devices.

This device was tested with accessories from the accessory list for electromagnetic compatibility. Other accessories can be used only if the electromagnetic compatibility is not interfered with. Use of noncompliant accessories can cause amplified electromagnetic emissions or decreased electromagnetic interference resistance in the device.

 WARNING
Hazard due to inadequate safety distance If high-frequency mobile communication devices are used to close to this device, malfunctions can occur that may endanger the patient. A safety distance of at least 0.3 m (1 foot) must be maintained.

Electromagnetic environment

The device is only to be used in environments indicated in the "Intended use" section of the operating manual. The medical device is intended only for operation in the electromagnetic environment indicated below.

Emissions	Correspond to	Electromagnetic environment
HF emissions EN 55011 (CISPR 11) Radiated: 30 MHz to 1 GHz Conducted: 150 kHz to 30 MHz	Class B, Group 1	The medical device is intended for use in all facilities, including residential buildings and facilities that are directly (without a transformer) connected to the same low voltage network as the residential building.
Harmonic emissions (IEC 61000-3-2)	Class A	
Emissions from Voltage fluctuations/flickering (IEC 61000-3-3)	Requirement is met.	

Immunity against	Test level and electromagnetic environment to be maintained	Electromagnetic environment
Electrostatic discharge (IEC 61000-4-2)	Contact discharge: ± 8 kV Air discharge: ± 15 kV	Floors made of wood, concrete or ceramic tile. For a synthetic floor covering, the relative humidity should be at least 30%.
Fast transient electrical disturbances/bursts (IEC 61000-4-4)	Power cable: ± 2 kV Longer signal input lines/signal output lines: ± 1 kV	The quality of the supply voltage should correspond to that of a typical business or hospital environment.
Impulse voltage/surges (IEC 61000-4-5)	Voltage: External conductor against external conductor: ± 1 kV External conductor against ground cable: ± 2 kV	
Voltage dips and short interruptions of the supply voltage (IEC 61000-4-11)	30% to 100%, 10 ms to 5 s, various phase angles	
Magnetic field at the supply frequency (IEC 61000-4-8)	50 Hz and 60 Hz: 30 A/m	Devices with unusually strong line-frequency magnetic fields (transformer stations, etc.) should not be operated near the medical device.
Emitted HF disturbance (IEC 61000-4-3)	80 MHz to 2.7 GHz: 10 V/m	Near equipment marked with the following symbol, disturbances are possible: 
Conducted HF interference (IEC 61000-4-6)	150 kHz to 80 MHz: $3 V_{rms}$ ISM and amateur radio bands: $6 V_{rms}$	

Recommended safety distances from portable and mobile HF communications equipment

Rated power of transmitter [W]	150 kHz – 800 MHz $d = 1.2 \sqrt{P}$	800 MHz – 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12 m (0.39 ft)	0.23 m (0.76 ft)
0.1	0.38 m (1.25 ft)	0.73 m (2.4 ft)
1	1.2 m (3.9 ft)	2.3 m (7.6 ft)
10	3.8 m (12.5 ft)	7.3 m (23.9 ft)
100	12m (39 ft)	23 m (76 ft)

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Waldmann Ljusteknik AB