Montage- und Gebrauchsanweisung	04	<b>DEU</b>
Installation and operating instructions	19	<b>ENG</b>
Mode d'emploi	35	<b>FRA</b>
Istruzioni per il montaggio e l'uso	50	<b>ITA</b>
Instrucciones de montaje y funcionamiento	66	<b>ESP</b>
Instruções de instalação e de utilização	81	<b>POR</b>
Montage- en gebruiksaanwijzing	97	<b>NLD</b>
Bruksanvisning	113	<b>SWE</b>
取り付けおよび使用の説明書	128	<b>JAP</b>

D <sup>med</sup> ®	<b>Triango</b>	100 C
D <sup>med</sup> ®	<b>Triango</b>	100 W
D <sup>med</sup> ®	<b>Triango</b>	100 F

Behandlungsleuchte  
 Treatment light  
 Lumière de traitement  
 Luce di trattamento  
 Luz de tratamiento  
 Candeeiro de Tratamento  
 Behandelingslicht  
 Manöverlampe  
 処置ライト

## SYMBOLE; SYMBOLS; SYMBOLES; SIMBOLI; SÍMBOLOS; SÍMBOLOS; SYMBOLEN; SYMBOLER; 記号の意味

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**IMPORTANT!**  
**THIS OPERATING MANUAL MUST BE THOROUGHLY READ BEFORE**  
**USING THE PRODUCT!**


→ **STORE FOR LATER USE!**

## CONTENT

1.	VERSIONS AND DELIVERY CONTENT	20
1.1	Triango 100 C.....	20
1.2	Triango 100 W.....	20
1.3	Triango 100 F.....	20
2.	SAFETY INSTRUCTIONS	21
2.1	Intended use.....	21
2.2	User profiles.....	21
2.3	Safety instructions.....	21
2.4	Warning levels.....	21
2.5	Specific fastening instructions.....	22
3.	INSTALLATION: Triango 100 C	22
3.1	Workload data.....	22
3.2	Shortening the ceiling tube.....	22
3.3.	Installing the ceiling holder.....	22
3.4.	Installing the ceiling tube.....	23
3.5	Installing the ceiling arm.....	24
4.	INSTALLATION: Triango 100 W	24
4.1	Workload data.....	24
4.2.	Installing the wall holder.....	24
4.3	Installing the wall mounting arm.....	25
5.	INSTALLATION: Triango 100 F	26
6.	INSTALLATION: luminaire head	27
7.	OPERATION	28
8.	CLEANING AND DISINFECTION	28
8.1	Sterilizing the hand grip.....	29
9.	SAFETY INSPECTIONS	29
9.1	SETTING THE SPRING FORCE.....	29
10	DISMANTLING	30
10.1	Disposal.....	30
11.	ACCESSORIES	30
12	ADDITIONAL INFORMATION	30
13.	TROUBLESHOOTING	31
14.	TECHNICAL DATA	31
15.	ELECTROMAGNETIC COMPATIBILITY (EMC)	33

# 1. VERSIONS AND DELIVERY CONTENT


## 1.1 Triango 100 C

**A:** Ceiling bracket, ceiling bonnet and end ring  
**C:** Lamp head with sterilisable handle

**B:** Ceiling pipe  
**D:** Ceiling arm


## 1.2 Triango 100 W

**A:** Wall bracket, wall bonnet and end ring  
**C:** Wall arm

**B:** Wall joint  
**D:** Lamp head with sterilisable handle

## 1.3 Triango 100 F

**A:** Lamp head with sterilisable handle  
**C:** Lower standpipe

**B:** Upper stand tube with spring arm  
**D:** Rolling frame

**E:** Mains cable

## 2. SAFETY INSTRUCTIONS

### 2.1 Intended use

The Triango 100 is a treatment luminaire (minor surgical luminaire). It is an individual luminaire for use in operating rooms that is used near the patient to support diagnosis or treatment, which poses no risk to the patient in the event of any interruption caused by light failure. It is meant for continuous operation and is not intended to be combined with other medical devices.

### 2.2 User profiles

#### Health Professional

All individuals who have completed medical training and work in the professional field they trained for.

#### Cleaning specialist

Trained in national and workplace hygiene regulations.

#### Qualified electrician

Trained in electronics and electrical technology and knows the relevant standards and regulations.

#### Qualified specialist

Qualified due to his technical training, knowledge and experience and knowledge of the rules, to perform assembly/dismantling.

### 2.3 Safety instructions

- ▶ Operation by health professionals
- ▶ This manual is part of the product and must be stored and made available to all future users.
- ▶ All work on the luminaire (including repairs) may only be performed by a qualified electrician. Assembly may only be done by a qualified specialist.
- ▶ The luminaire may not be altered or tampered with. Only approved original parts may be used. Any use other than that intended, using original parts, can change technical parameters and pose a death hazard.
- ▶ Do not exceed the maximum weight and do not hang, lean or climb on the arm, because this may cause the device to tip over, which could result in serious injury.
- ▶ Operation in potentially explosive areas is prohibited. The luminaire's current supply is a potential ignition source.
- ▶ The luminaire must be used only in dry, dust-free rooms.
- ▶ The luminaire should not be left on without supervision.
- ▶ Only connect the luminaire to the supply network with a protective ground conductor (PE), to prevent electrical shock.
- ▶ In the case of luminaires of protection class I, the protective ground conductor (PE) must be connected to luminaires housing.
- ▶ Do not use a damaged luminaire. Defective cables and a defective hand grip also pose a potential hazard. Do not place a cord near heat sources or sharp edges.
- ▶ Do not place extra loads on the luminaire head and arm system.
- ▶ The luminaire must not be covered with a cloth or similar item during operation.
- ▶ The ventilation openings (if present) must always be kept clear during operation!
- ▶ The luminaire must not be operated near external heat sources that exceed the luminaire's maximum ambient temperature.
- ▶ The luminaire must not be used outside the specified ambient conditions.
- ▶ Do not use with medical devices that may react sensitively to a light spectrum within the visible range (such as pulsating light and/or light with high illumination intensity).
- ▶ The luminaire may only be used for the intended use described here.
- ▶ The manufacturer cannot be held responsible for any damages resulting from use deviating from the intended use, or failure to observe the safety instructions and warnings.
- ▶ When using more than one luminaire at the same time, the total illumination intensity must not exceed  $E_e 1000\text{W}/\text{m}^2$ .
- ▶ Before connecting to the power grid, it is essential to check that the grid data match the device data.

### 2.4 Warning levels

#### DANGER

Warning of hazards that can result in **death or serious injury** if instructions are not followed.

#### WARNING

Warning of hazards that can result in **injury** if instructions are not followed.

#### CAUTION

Warning of hazards that can cause **material damage** if instructions are not followed.

## 2.5 Specific fastening instructions

### Triango 100 W

- ▶ **Fastening material is not included with delivery.**
- ▶ The ceiling mounting must only be mounted on ceilings having a concrete stability class B25 (C20/25) or higher.
- ▶ Solid ceiling reinforcement parts must not come into contact with the ceiling mounting. If there is any doubt, an authorized professional must confirm that mounting is possible on the mounting base being used. The load-bearing capacity of the ceiling construction must be planned, checked and confirmed by a structural engineer beforehand.
- ▶ Borings must be made by specialists in compliance with the boring tolerances for the reinforcement anchor that have been approved by the manufacturer. In the event of faulty drilling, such as drilling through a metal support bar, a structural engineer must be called in.
- ▶ The luminaire should be installed in such a way as to avoid permanent strain on the height stops during operation.
- ▶ If there is plastering or cladding over the concrete, make sure the fixing anchor is properly attached right into the concrete.
- ▶ Screws must be tightened carefully using a torque wrench that complies with the instructions of the mounting element manufacturer.

### Triango 100 W


- ▶ **Fastening material is not included with delivery.**  
The luminaire must only be mounted on walls that can guarantee a secure hold. Specialist staff will find details of the requirements in Chapter 5.1 (Load data).

## 3. INSTALLATION: Triango 100 C


### 3.1 Workload data

Bending moment $M_B$	245 Nm
Vertical weight $F_G$	160 N

### 3.2 Shortening the ceiling tube

- ▶ Before shortening the ceiling tube, remove the cable from the ceiling tube.
- ▶ Use a hack saw to shorten the upper end of the ceiling tube to the desired length, and then deburr it.

- ▶ Remove the fastening screw "C".
- ▶ Insert the ceiling tube "A" into the ceiling bracket "B" and drill the ceiling holder's existing hole to 9 mm in diameter. Drill the opposite hole separately.
- ▶ Note: Pull the cable through from the lower side to the upper pipe side (the 3-pin plug first) after sawing and drilling.

### 3.3. Installing the ceiling holder

#### DANGER

##### Assembly by qualified personnel

- ▶ Assembly must be done by qualified personnel only. Lack of appropriate knowledge may pose a death hazard.
- ▶ Two people are needed for assembly.

#### DANGER


##### Death hazard from falling luminaire.

- ▶ The ceiling mount must only be mounted on ceilings with a concrete stability class of B25 (C20/25) or higher.
- ▶ There should be no contact with reinforcement parts of the solid ceiling during mounting. When in doubt, a licensed technician must approve installation on the specific installation surface. The load bearing capacity of the ceiling structure must be planned, checked and confirmed by a structural engineer beforehand.
- ▶ The holes must be drilled by a professional and in compliance with the drilling tolerances permitted by the manufacturer of the fastening anchor. If a drilling error occurs, such as drilling into a reinforcement bar, a structural engineer must be consulted.
- ▶ Install the luminaire such that the vertical stops are not continuously stressed during operation.
- ▶ With plaster or cladding in front of concrete, the fastening anchor must be driven completely into the concrete.
- ▶ The screws must be tightened carefully using a torque wrench according to the specifications of the fastening anchor manufacturer.

#### DANGER

##### Death hazard from electric shock.


- ▶ All poles of the power cable must be disconnected from the power grid from an external lockable switch.


- ▶ Draw 6 drill marks.
- ▶ Observe the position of opening "A" for connecting the current.

**⚠ WARNING**


Wear safety equipment according to the tool manufacturer's instructions.

- ▶ Drill holes and blow out with a bellows.

- ▶ Check the hole distances.


- ▶ Hold the ceiling bracket to the ceiling and hammer in the fastening anchor.
- ▶ Tighten the fastener according to manufacturer's instructions.

**⚠ DANGER**


**Death hazard from electric shock.**

- ▶ Do not turn on the power supply before the luminaire head has been installed.
- ▶ This device must only be connected to a power grid with a protective ground conductor (PE) to prevent risk of electric shock


- ▶ Connect the power.


**3.4. Installing the ceiling tube**


- ▶ Pull the ceiling tube cable through the opening "B" of the ceiling holder.
- ▶ Insert the ceiling tube "A" into the ceiling holder.
- ▶ Secure the safety screw "C" and the M8 nut.

- ▶ Tighten the M8 safety screw "C" and the nut (20 Nm).
- ▶ Tighten all 4 grub screws "D" (5 Nm).

- ▶ Connect the ceiling tube plug with the power supply plug.


- ▶ Slide the ceiling cover and ring over the ceiling bracket and screw them tight.

### 3.5 Installing the ceiling arm

**WARNING**


The opening spring balanced arm is an injury hazard

- ▶ Do not remove the transport safeguard until after the luminaire head has been mounted


- ▶ Connect the plug from the ceiling tube with plug on the ceiling arm
- ▶ Carefully push protruding cable into the ceiling tube

Perform this assembly step on an even surface:

- ▶ Insert the sliding piece 1 above into the ceiling arm groove 2
- ▶ Push the ceiling tube 3 and ceiling arm 2 together
- ▶ Secure the sliding piece 1 by screwing in the threaded stud 4 (up to the stop)
- ▶ Loosen the threaded stud 4 again by a maximum of a ¼ turn and checking the rotating function.

- ▶ Continue with luminaire head assembly, see Chap. 6

## 4. INSTALLATION: Triango 100 W

### 4.1 Workload data

Bending moment $M_B$	275 Nm
Vertical weight $F_G$	155 N

### 4.2. Installing the wall holder

- ▶ Fastening hardware does not come with delivery.

**DANGER**

Assembly by qualified personnel

- ▶ Assembly must be done by qualified personnel only. Lack of appropriate knowledge may pose a death hazard.
- ▶ Two people are needed for assembly.

**CAUTION**

Determine the fasteners to use from the workload data table

- ▶ Before assembly, make sure the rod dimensions are correct.

**CAUTION**


Make sure the wall holder is in the correct position.

- ▶ The wall holder must be aligned according to the axis in the image.
- ▶ Misalignment will compromise mechanical safety.
- ▶ We recommend using a counterplate on lightweight construction walls (not included with delivery).

**DANGER**

Death hazard from electric shock.


- ▶ All poles of the luminaire must be disconnected from the power supply by an external lockable switch.


- ▶ Draw 6 drill marks.
- ▶ Observe the position of opening "A" for connecting the current.

**WARNING**


Wear safety equipment according to the tool manufacturer's instructions.

▶ Drill holes and blow out with a bellows.

▶ Check the hole distances.


- ▶ Position the wall bracket on the wall and hammer in the fastening anchor.
- ▶ Tighten the fastener according to manufacturer's instructions.


**⚠ DANGER**

**Death hazard from electric shock.**


- ▶ Do not turn on the power supply before the luminaire head has been installed.
- ▶ This device must only be connected to a power grid with a protective ground conductor (PE) to prevent risk of electric shock


▶ Connect the power.


▶ Insert wall bracket "A" (with end ring "C" and cover "B" in position) into the wall holder "D" and simultaneously pull the plug through the rectangular cutout.

- ▶ When the wall bracket is vertically aligned, install it with the M8 locking screw and nut "E" and tighten (20 Nm).
- ▶ Tighten all 4 grub screws "F" (5 Nm).

▶ Connect the ceiling tube plug with the power supply plug.


▶ Push the cover against the wall and screw tight with the ring (0.5 Nm)

**4.3 Installing the wall mounting arm**


**⚠ WARNING**

**The opening spring balanced arm is an injury hazard**


- ▶ Do not remove the transport safeguard until after the luminaire head has been mounted.


▶ Connect the wall-mounted joint plug to the wall bracket plug


▶ Connect the wall arm plug to the wall-mounted joint plug

- ▶ Insert the arm for wall-mounting **3** into wall bracket **1** slightly
- ▶ Insert the sliding piece **2** into the groove
- ▶ Push up the arm **3** in such a way that the hole in the wall bracket is aligned with the hole in the sliding piece
- ▶ Screw grub screw **4** into wall bracket


- ▶ In the case of a dead stop, screw this back by a ¼ turn and check the turning function


- ▶ Continue with luminaire head assembly, see Chap.6


## 5. INSTALLATION: Triango 100 F


- ▶ Remove the M3 screws on the sides and remove the power supply tray "A".


- ▶ Loosen the 2 screws on the lower support tube.
- ▶ Feed the support tube cable through the roller stand.
- ▶ Attach the support tube to the roller stand using the 2 Allen screws and serrated washers (10 Nm).
- ▶ The cable holder must be aligned rearward.


- ▶ Connect the same-colored stranded wires of the lower support tube and of the power supply.

- ▶ Reinstall the power supply tray with the M3 Allen screws and serrated washers.


- ▶ Loosen the grub screw **1**
- ▶ Place the spring balanced arm **2** on the upper support tube
- ▶ Attention: Grub screw **1** must be inserted into slot **3**.


- ▶ Screw in the grub screw in the alignment shown
- ▶ In the case of a grub screw dead stop, screw this back by a ¼ turn and check the turning function (Swing range ± 30°)

- ▶ Push together the plugs on the upper and lower support tube until they click.
- ▶ Place the upper support tube on the lower one.

- ▶ Tighten the cable holder using a type 3 Allen screw on the upper support tube (2.4 Nm).

**⚠ WARNING**

**Never unscrew the upper cable holder — injury hazard.**

- ▶ If both cable holders are unscrewed, the connecting piece will come loose and fall down, which can cause personal injury and damage to the cable and device.

**⚠ DANGER**

**Death hazard from electric shock.**

- ▶ Do not turn on the power supply before the luminaire head has been installed.

**6. INSTALLATION: luminaire head**

**⚠ DANGER**

**Death hazard from electric shock.**


- ▶ Before any work on the device, disconnect the current supply or pull the plug out of the socket and secure the device against being turned back on.

**⚠ WARNING**


**Injury hazard**

- ▶ The spring balanced arm is under high spring stress. If the boom is lowered without the installed device, it **MUST** be held firmly.
- ▶ If it is let go, it springs upward and can cause severe injury.
- ▶ Never remove the luminaire body unless the boom is in the upper position or is held securely in low position by a second person.


- ▶ To prevent severe injury or damage, always get the help of a second person when installing and removing the luminaire body.

- ▶ Connect the luminaire head plug to the wall-mounted joint plug

- ▶ Before assembly one must ensure that the brass nose (see arrow) and the cardan joint bracket are on the left


- ▶ Push together the luminaire head 1 and arm 3; in doing so the cardan joint bracket must be on the left 2 of the luminaire head
- ▶ then screw in both countersunk screws M4 (both with a toothed washer) in (1.5Nm)
- ▶ Check the turning function

**⚠ WARNING**

**Injury hazard from falling luminaire head.**

- ▶ Make sure the locking device is properly installed.
- ▶ Damage to due to a poorly mounted luminaire head.
- ▶ Tighten the screw as the luminaire head's friction requires.


- ▶ Attach the hand grip.

## 7. OPERATION

### CAUTION

- ▶ This product emits potentially hazardous optical radiation. Do not stare directly into the light cone. Eye irritation may occur.
- ▶ The radiation emitted by this product complies with the exposure limits for reducing the risk of photobiological hazards based on IEC 62471: RG 2 (medium risk).

- ▶ Switch the luminaire on or off using the button 1.
- ▶ The luminaire can be dimmed using the +/- buttons 2
- ▶ The light color can be adjusted using button 3
- ▶ The focus can be adjusted using button 4
- ▶ The endoscope mode can be adjusted using button 5
- ▶ Before each use, perform a function test: all LEDs in the light cone must come on.

### Triango 100 F

#### DANGER

##### Electric shocks are life-threatening.

- ▶ Do not plug in any damaged power cables.
- ▶ If there are any signs of damage to the power cable, immediately replace it with a new one.
- ▶ Connection voltage and frequency must match data on the type plate.
- ▶ Only connect the luminaire to the supply network with a protective earth conductor (PE).

#### DANGER


##### In the event of a change of location

- ▶ Put the terminal into its lowest position.
- ▶ Unlock the castors.
- ▶ Take care not to run over anything, including the connection cable.
- ▶ Hold firmly during transport within the hospital.
- ▶ Take particular care on slopes, thresholds, bumps and other obstacles.

### CAUTION

#### Triango 100 F

- ▶ Wind the power cable onto the cable holder when not in operation

- ▶ Plug in cable
- ▶ Connect cable to the network
- ▶ Lock the castors

## 8. CLEANING AND DISINFECTION

### DANGER

#### Death hazard from electric shock

- ▶ Before disinfection cleaning, pull the plug out of the socket and secure the device against being turned back on.

### CAUTION

#### Material damage due to incorrect cleaning

- ▶ For cleaning, use only products that do not impair operation of the luminaire.
- ▶ For cleaning, do not use solvent-based, chlorine-based or abrasive detergents, because they can crack plastic parts and cause other damage.
- ▶ The cleaning agents must be approved for use on plastics such as PC, PMMA, PA and ABS.
- ▶ Concentrated disinfectant can damage the luminaire.
- ▶ For concentration and application times, check the information provided with the product used.
- ▶ The wrong cloth may cause scratches.


### RECOMMENDED DISINFECTANTS

- ▶ Bacillol 30 Foam
- ▶ Dismozon Plus
- ▶ Kohrsolin Extra
- ▶ Lysoformin
- ▶ Mediclean
- ▶ Meliseptol Foam Pure
- ▶ Microbac Tissues
- ▶ Mikrozid Sensitive Liquid
- ▶ Mikrozid AF Liquid
- ▶ Sagrotan
- ▶ Terralin Protect
- ▶ Virex Tb
- ▶ CaviCide 1
- ▶ Germicidal Bleach
- ▶ Hexaquart XL

**CAUTION**


**Dirt reduces the light strength**

- ▶ Keep cover clear through regular cleaning.
- ▶ Only wipe cleaning allowed.

- ▶ Clean the PA clear cover using a non-abrasive cleaning cloth and a suitable cleaning agent

**CAUTION**

- ▶ Wipe cleaning is only allowed in a horizontal position.

**CAUTION**

To minimize the risk of disease transmission, applicable health and safety regulations and the requirements of the national hygiene and disinfection authorities must be observed in addition to these instructions.

**8.1 Sterilizing the hand grip**

- ▶ Sterilization must be done according to **ISO 17665-1** (Sterilization of health care products in moist heat).

**CAUTION**

**Damage to the hand grip**

- ▶ Do not sterilize with hot air
- ▶ Package the hand grip in a sterile bag before sterilization.
- ▶ The hand grip is designed exclusively for damp sterilization with 3 times fractioned pre-vacuum and saturated steam with the following parameters:

Temperature	134°C
Overpressure	2.0 bar
Dwell time	6 min
Vacuum drying	20 min

- ▶ After sterilization, check the hand grip for mechanical integrity.
- ▶ Do not use damaged hand grips.

**9. SAFETY INSPECTIONS**

**⚠ DANGER**

**Death hazard from electric shock.**

- ▶ Unplug from the grid and turn switch to off position.
- ▶ The power supply cable must be checked at least once a year for damage.

**CAUTION**

- ▶ Maintenance and repairs can be performed only by qualified electricians.
- ▶ The corresponding user profile is in Chapter 2 Safety instructions.


**EVERY YEAR:**

- ▶ Check the power supply cable for damage and replace if necessary.
- ▶ Check for deformations and cracks in the plastic parts.
- ▶ Check the load bearing system for distortion or damage.
- ▶ Check for loose parts.

**9.1 SETTING THE SPRING FORCE**

**CAUTION**


- ▶ The ideal spring force is set at the factory.

- ▶ As a first measure slightly tighten the joint cover (by a ¼ turn); this to achieve greater stability of the luminaire head, if necessary

Otherwise:

- ▶ Remove the service cover of the spring balanced arm (loosen the 2 crosstip screws)

- ▶ Place a suitable tool (e.g. a 3-type Allen key) into the hole 1, turn by ¼ turn in the direction of the arrow, a repeat a number of times until the desired stability is reached
- ▶ Screw the service cover on again (attention: do not tighten too strongly)

## 10 DISMANTLING

**⚠ DANGER**

**Death hazard from electric shock.**

Before dismantling the luminaire, disconnect all pins from the power supply.

**⚠ WARNING**

**Injury hazard**

The spring balanced arm is under a high spring force. If the terminal device is not in the uppermost spring balanced arm position, the spring balanced arm will move rapidly upward and cause severe injury. Only dismantle the terminal device when the spring balanced arm is in the uppermost position

### 10.1 Disposal

Do not discard the luminaire with household waste. Follow local regulations and bring the luminaire to a disposal site or give it to a dealer with appropriate service.  
Cut the cable directly at the casing.

 The products listed above are over 95% recyclable. For a high percentage of the used materials to be either physically reused or used for energy after the end of this product's life cycle, the luminaires have been designed with recycling in mind. They do not contain materials that are hazardous or require monitoring.

## 11. ACCESSORIES

► Ergonomic sterilizable handle (order no D10.442.000)

► Ergonomic sterilizable handle (order no D10.442.000)

## 12 ADDITIONAL INFORMATION

The luminaire itself is maintenance free.

Additional documents may be requested from the manufacturer for this product.

Using this luminaire does not pose a risk to other equipment. To save energy, the luminaire should be switched on only when it is actually needed.

Any serious incident that has occurred with the product **must be reported** to the manufacturer or their representative and the responsible authorities of the country in which the user is located.

### 13. TROUBLESHOOTING

Fault	Possible cause	Troubleshooting	User profile
The luminaire does not light up	Contact fault	Switch on again	All
The luminaire does not light up	Lamp defective	Contact manufacturer's service dept.	Only by manufacturer's service dept.
The luminaire does not light up	No mains voltage	Check voltage, check all connections	Qualified Electrician

### 14. TECHNICAL DATA

Electrical Data:	
Rated input voltage	100 – 240 VAC
Frequency range	50 / 60 Hz
Power consumption	61 – 80 VA
Input current	610 – 330 mA
Power factor	0.69 – 0.45
Integral electronic transformer	24 VDC output
Photometric values:	
Central illuminance $E_c$ at 1.0m distance	100'000 lx **
Light field diameter $d_{10}$ at a distance of 1.0 m: Triango 100-1, Triango 100-3, Triango Endo 100-1 Triango Fokus 100-1, Triango Fokus 100-3	dia. = 18 cm * dia. = 18 cm / 23 cm / 28 cm*
Light field diameter $d_{50}$ at a distance of 1.0 m:	dia. = 10 cm *
Color temperature: Triango 100-1, Triango Fokus 100-1, Triango Endo 100-1 Triango 100-3, Triango Fokus 100-3	4300K * 3700K*/ 4300K*/ 4700K*
Color rendering Index $R_a$	$\geq 96^*$
Color rendering Index $R_9$ : Triango 100-1, Triango Fokus 100-1, Triango Endo 100-1 Triango 100-3, Triango Fokus 100-3	90* 95*
Depth of illumination $L_1 + L_2$	80.4 cm *
Total irradiance $E_e$ :	415 W/m <sup>2</sup>
Relationship between irradiance $E_e$ and irradiance $E_c$ :	3.51 mW/m <sup>2</sup> /lx
Remaining irradiance:	1 shade: 27.9% 2 shades 0°: 63.4% 2 shades 45°: 48.5% 2 shades 90°: 25.3% 2 shades 135°: 62.3% Tube: 100% Tube + 1 shade: 28.2% Tube + 2 shades 0°: 62.6% Tube + 2 shades 45°: 49.1% Tube + 2 shades 90°: 25.2% Tube + 2 shades 135°: 60.5%
	* -10% / +20% tolerance
	** only for 4300K / $D_{10} = 180\text{mm}$

<b>Ambient conditions for transport, storage and operation:</b>	
Ambient temperature (storage and transport)	-20°C to +70°C
Ambient temperature (operation)	10°C to +35°C
rel. air humidity (non-condensing)	max. 75%
<b>Weight:</b>	
Luminaire head	2.4 kg
<b>Triango 100 C</b>	15.9 kg
<b>Triango 100 W</b>	15.3 kg
<b>Triango 100 F</b>	20.0 kg
<b>Operating mode:</b>	
Operating mode	Continuous operation
<b>Classification:</b>	
Triango 100	Protection class I
Degree of protection according to IEC 60529	IP 20
Luminaire head	IP 43 (horizontal position)
Classification according to EU-REGULATION 2017/745 (MDR), Article 51	Class I
U.S. FDA Device Class	Class I
Electrical safety testing and EMC according to:	AAMI ES60601-1 : 2005/A2 :2010/(R)2012 CAN/CSA-C22.2 No 60601-1:14 IEC 60601-1:2005 + A1 :2012 IEC 60601-1-2:2014 IEC 60601-2-41:200- + A1 :2013
Blue light danger according to EN/IEC 62471	RG 2 (moderate risk)
<b>Life cycle of the light source:</b>	
Life cycle	50,000 h (L70/B50)

## 15. ELECTROMAGNETIC COMPATIBILITY (EMC)

Electrical medical devices are subject to special precautionary measures regarding electromagnetic compatibility. This device can be affected by other electrical devices.

This device was tested with accessories from the accessory list for electromagnetic compatibility. Other accessories can be used only if the electromagnetic compatibility is not interfered with. Use of noncompliant accessories can cause amplified electromagnetic emissions or decreased electromagnetic interference resistance in the device.

### WARNING

#### Hazard due to inadequate safety distance

If high-frequency mobile communication devices are used too close to this device, malfunctions can occur that may endanger the patient. A safety distance of at least 0.3 m (1 foot) must be maintained.

### Electromagnetic environment

The device is only to be used in environments indicated in the "Intended use" section of the operating manual.

The medical device is intended only for operation in the electromagnetic environment indicated below.

Emissions	Correspond to	Electromagnetic environment
HF emissions EN 55011 (CISPR 11) Radiated: 30 MHz to 1 GHz Conducted: 150 kHz to 30 MHz	Class B, Group 1	The medical device is intended for use in all facilities, including residential buildings and facilities that are directly (without a transformer) connected to the same low voltage network as the residential building.
Harmonic emissions (IEC 61000-3-2)	Class A	
Emissions from Voltage fluctuations/flickering (IEC 61000-3-3)	Requirement is met.	

Immunity against	Test level and electromagnetic environment to be maintained	Electromagnetic environment
Electrostatic discharge (IEC 61000-4-2)	Contact discharge: $\pm 8$ kV Air discharge: $\pm 15$ kV	Floors made of wood, concrete or ceramic tile. For a synthetic floor covering, the relative humidity should be at least 30%.
Fast transient electrical disturbances/bursts (IEC 61000-4-4)	Power cable: $\pm 2$ kV Longer signal input lines/signal output lines: $\pm 1$ kV	The quality of the supply voltage should correspond to that of a typical business or hospital environment.
Impulse voltage/surges (IEC 61000-4-5)	Voltage: External conductor against external conductor: $\pm 1$ kV External conductor against ground cable: $\pm 2$ kV	
Voltage dips and short interruptions of the supply voltage (IEC 61000-4-11)	30% to 100%, 10 ms to 5 s, various phase angles	
Magnetic field at the supply frequency (IEC 61000-4-8)	50 Hz and 60 Hz: 30 A/m	Devices with unusually strong line-frequency magnetic fields (transformer stations, etc.) should not be operated near the medical device.
Emitted HF disturbance (IEC 61000-4-3)	80 MHz to 2.7 GHz: 10 V/m	Near equipment marked with the following symbol, disturbances are possible: 
Conducted HF interference (IEC 61000-4-6)	150 kHz to 80 MHz: 3 V <sub>rms</sub> ISM and amateur radio bands: 6 V <sub>rms</sub>	

<b>Recommended safety distances from portable and mobile HF communications equipment</b>		
<b>Rated power of transmitter [W]</b>	<b>150 kHz – 800 MHz <math>d = 1.2 \sqrt{p}</math></b>	<b>800 MHz – 2.5 GHz <math>d = 2.3 \sqrt{p}</math></b>
0.01	0.12 m (0.39 ft)	0.23 m (0.76 ft)
0.1	0.38 m (1.25 ft)	0.73 m (2.4 ft)
1	1.2 m (3.9 ft)	2.3 m (7.6 ft)
10	3.8 m (12.5 ft)	7.3 m (23.9 ft)
100	12m (39 ft)	23 m (76 ft)

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Waldmann Eclairage S.A.S
### Derungs Licht AG
### Derungs Medical GmbH